

Hemodynamic Monitoring: The Septic Patient

Clinical Profile:

Sepsis (Sĕp'sis) [Gr. sĕpsis, putrefaction]. Pathologic state, usually febrile, resulting from the presence of microorganisms or their poisonous products in the bloodstream.¹

Sepsis must be assessed, treated and arrested utilizing accurate and efficient diagnostic modalities.² Given the marked right ventricular (**RV**) dysfunction demonstrated by septic and hypovolemic patients, preload measurement most accurately was assessed, historically, by using the filling pressures, central venous pressure (**CVP**) and pulmonary artery wedge pressure (**PAWP**)⁷, due to an assumed correlation with diastolic volume.

Correlations:

- **Correlation #1: between PCWP and LVEDV**

It has been shown that pulmonary capillary wedge pressure (**PCWP**) did not accurately reflect changes in left ventricular (**LV**) preload. Left ventricular end diastolic volume (**LVEDV**) was found to depend on right ventricular (**RV**) function and pulmonary vascular resistance (**PVR**). In patients with a positive response to volume loading, there was an increase in either **LVEDV** or left ventricular ejection fraction (**LVEF**).³

- **Correlation #2: between RVEDVI and CI**

Right ventricular end diastolic volume index (**RVEDVI**) has been shown to be a better predictor of cardiac index (**CI**) ($r=0.61$) than pulmonary artery wedge pressure (**PAWP**) ($r=0.42$). Discrepancies between **PAWP** and **RVEDVI** occurred in more than 50% of patients studied. “**RVEDVI** more accurately predicted preload recruitable increases in cardiac output (**CO**).”⁴

- **Effect of monitoring in cases of severe septic shock**

It has been shown that “**RV** monitoring during severe septic shock permits the identification of patients with **RV** dysfunction who are difficult to identify with usual measurements.”⁵

- **Correlation #3: between RVEF interpretation and indication for treatment of septic patients**

Right ventricular ejection fraction (**RVEF**) values should

be interpreted after accounting for the increase in **RV** afterload and the inverse relationship between **RVEF** and pulmonary artery (**PA**) pressures. In acute respiratory failure, severe sepsis or trauma and during or after surgery, these measurements have been used to assess **RV** function.

Despite a well preserved stroke volume (**SV**) and **CI**, septic shock patients can have a low **RVEF**. A recent series of these patients that were studied by Jean-Louis Vincent, M.D. demonstrated a reduction in **RVEF** early in the course of septic shock with a return to baseline **RVEF** during the resolution of shock. Dr. Vincent concluded that, prognostically, a very low **RVEF** was a bad sign. Vasopressors may be indicated for treatment in septic patients when arterial hypotension and pulmonary hypertension threaten **RV** perfusion.⁶

Outcome: Effect of combining SvO₂, REF, EDV and CCO in a single catheter

It has been recognized that a new catheter, combining continuous venous oximetry, **RVEF**, end diastolic volume (**EDV**) measurements, and continuous thermodilution cardiac output (**CTCO**), will yield not only information regarding oxygen consumption and delivery, but will also allow an even better understanding of hemodynamics. A continuous indication of pulmonary vascular resistance (**PVR**) and systemic vascular resistance (**SVR**) may be used as an index to left ventricular (**LV**) and right ventricular (**RV**) afterload when continuous cardiac output (**CCO**) is integrated with continuous central and arterial pressure measurements.⁷

Loren D. Nelson, M.D. concluded in his report that the combination of intravascular pressure measurements and **CCO** also allows calculation of **LV** and **RV** stroke work, “that, under conditions of steady state preload and afterload, reflects the contractile function of the heart as a determinant of cardiac performance.”⁷ On a near continuous basis, the critical care team will thus be provided with complete hemodynamics and oxygen transport values.



Edwards
LIFESCIENCES

The quotes and information referenced in this brochure were selected from independent third party publications and are not intended to suggest that such third parties have reviewed and/or endorsed Edwards' products.

References:

1. Taber, CW (1977) Taber's Cyclopedic Medical Dictionary, (13th edition) Philadelphia: FA Davis Company.
2. Price SA and Wilson LM, "Pathophysiology, Clinical Concepts of Disease Processes", 1982; 6:64-71
3. Hoffman MJ, et al., "Unsuspected Right Ventricular Dysfunction in Shock and Sepsis", Ann Surg, 1983; 198:307-319.
4. Diebel LN, et al., "End-Diastolic Volume: A Better Indicator of Preload in the Critically Ill", Arch Surg, 1992; 127: 817-821; Discussion 821-822.
5. Le Tulzo Y, et al., "Effects of Epinephrine on Right Ventricular Function in Patients with Severe Septic Shock and Right Ventricular Failure: A Preliminary Descriptive Study", Intensive Care Med, 1997; 23:664-670.
6. Vincent JL, "The Measurement of Right Ventricular Ejection Fraction", Intensive Care World, 1990; 7:133-136.
7. Nelson LD, "The New Pulmonary Arterial Catheters: Right Ventricular Ejection Fraction and Continuous Cardiac Output", Critical Care Clinics, 1996; 12:795-818.



Certified according
to European MDD

See instructions for use for full prescribing information.

Edwards Lifesciences, Edwards, and the stylized E logo are trademarks of Edwards Lifesciences Corporation.

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

© Copyright 2002 Edwards Lifesciences LLC
All rights reserved. 3164-12/01-CC



Edwards Lifesciences